



Analytical Technologies Birthing Table-



Specifications:-

ATL 5200 BT

A. Description of Function

1. Delivery bed is used for Baby Delivery and incorporate ideal blend of the patient's individual requirements on comfort and the professional needs of the delivery team, focusing on the aesthetic and functional design of the entire product.

B. Operational Requirements

1. Delivery bed is supplied with all accessories

C. Delivery Bed should have following essential specifications:

1. It have control devise for making height (44cm to 90cm) and back adjustments.[manual as well as remote control].
2. It have collapsible side rails
3. It have three sectional mattress and seat section should have large perineal cut. The mattress thickness is 50mm or more.
4. Head board and food section can be detached or slides and stores under the bed.
5. Wheels (dia- 6" or 8") provided with locking system.
6. Retractable foot section with indication for locking, so as to convert bed into table.
7. Infusion rods which have adjustable heights, quick release and attaches to all corners of bed.
8. Adjustable leg rests available as an accessory
9. Push grip handles
10. Sliding stainless steel bowl at perineal part of table
11. Catheter bag holder which can be attached on either side of bed

12. Able to give trendelenburg, reverse trendelburg and 60 degree sitting position both mechanically and electronically.
13. Adjustable foot supports for nursing staff
14. Easy to clean, sterilize (especially blood stains) and maintain
15. Frame is of epoxy powder coated steel
16. Dimensions - Length: Minimum 180 cm and width: Minimum 75 cm
17. Pelvic tilt: – 15 degree.
18. Easy slide calf supports swing into correct positional lock with single lever.
19. CPR release.
20. Weight capacity: 200 (Approx)



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395
Email: info@ais-india.com, info@hplctechnologies.com, info@analyticalgroup.net
Web: www.ais-india.com / www.analyticalgroup.net / www.hplctechnologies.com

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement

Analytical Technologies Limited.

www.analyticalgroup.net