

Blood Weighing Scale With Monitor-BWS 3000 Series



SINGLE PAN SCALE

Special Features:

LCD display, displays the weight and volume with an accuracy of 1 gm/ml.

Helps better balancing of refrigerated centrifuge.

Easy conversion of weight to volume.

Tare provision to account for the weight of the blood bag

Compact model BWS molded body & Auto calibration

Technical Specification:

Catalogue No.	BWS-3000 Series
Input voltage	230 VAC
Power supply	SMPS
Transducer	Load Cell
Weighing Range	3 Kg
Accuracy	1gm/ml
Calibration	Auto
Display	LCD with back light
Weight	2 kg with battery
Dimension	175 x 225 mm
Pan	SS 304
Body	ABS Mould
Vision	Audio & Visual Alarm
Volume	Set Volume