



Analytical Technologies Bronchoscope



BRONCHO 3700 SERIES

Specifications:-

A. PAEDIATRIC BRONCHOSCOPES EXTENDED LENGTH WITH PROXIMAL ILLUMINATION Special features:

1. Adjustable proximal prismatic light deflector
2. Maximum lumen and small external diameter
3. No loss of lumen through light carrier
4. Pediatric bronchoscope –extended length
5. Consisting of: Bronchoscope sheaths of size 2.5Fr, 3Fr, 3.5 Fr, 4Fr, 4.5Fr and 5 Fr of standard length.
6. Prismatic light deflector with connection for fibre optic light cable.
7. Glass window plug
8. Rubber telescope guide
9. FLUVOG adapter with sliding glass window plug, sealing cap, notched lens and Keyhole opening, movable.
10. Injection cannula , O.D 3.5mm,
11. For positive pressure assisted ventilation system, for use with bronchoscopes
12. Instrument guide for aspiration catheter.
13. Adapter to respirator
14. Sealing plug for respiration connector.
15. Adjustable magnifier, swing –away type

B. ACCESSORIES:

1. Straightforward telescope 0°, diameter 2.8mm, length 44cms. Fibre optic light transmission incorporated. Colour code: green
2. Lateral telescope 30°, diameter 2.8mm,length 44 cm fibre optic light transmission. Colour code: yellow.
3. Pediatric Optical forceps for peanuts and soft foreign bodies, for use with mentioned telescopes.
4. Pediatric Optical alligator forceps, for use with telescopes.
5. Rigid suction tube with rubber tip, working length 50cm, straight and curved Sponge holder, spring holder, working length 50 cm,
6. Radio opaque plastic tubing for aspiration, 7 Fr, 5 m coil, with adapter.
7. Telescopic Bridge for a fixed position between telescope and bronchoscope.
8. Grasping Forceps with peanut, Crocodile/ alligator Jaw, single action, 1.5 mm shaft diameter, 45 cm length (pediatric)
9. Cup Grasping Forceps, 1.5 mm shaft diameter, 35 cm length
10. 1Forceps peanuts and soft foreign bodies, double action jaws, sheath diameter 1.5 mm working length 35 cm
11. Micro Scissors, 30 cm length, 3.0 mm shaft diameter (pediatric).



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395
Email: info@ais-india.com, info@hplctechnologies.com, info@analyticalgroup.net
Web: www.ais-india.com / www.analyticalgroup.net / www.hplctechnologies.com

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement

Analytical Technologies Limited.

www.analyticalgroup.net