



Analytical Technologies Colposcope



COLON 3380 SERIES

Specification:-

- A. Qualify stereoscopic optics. Facilities assessment of the finest epithelial changes.
- B. Ergonomic design: For convenient and fast positioning and focusing of colposcope
- C. Compact size: To conveniently move from one room to another
- D. Inclined Binocular tubes 0-180 degree or more: for best viewing experience
- E. Objective lenses: with different focal lengths allow user to select the convenient working distance. 300 mm, 400 mm
- F. Stepless zoom magnification settings: allow user to study epithelium at high magnification and carry you treatment at low magnification.
- G. Swing in vessel delineation filter: For improved visual activity.
- H. Optimum cold light illumination system: Help distinguishing small color differences in epithelium
- I. Suitable for tuboplasty (with swivel Arm Stand) & other Micro Surgical Procedure.



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395
Email: info@ais-india.com, info@hplctechnologies.com, info@analyticalgroup.net
Web: www.ais-india.com / www.analyticalgroup.net / www.hplctechnologies.com

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement