



Analytical Technologies DENTAL CHAIR (BASIC)



Specification:-

ATL 3030 DC

A. Description of Function Dental Chair is the dental chair required for dental examination and dental procedures.

B. Operational Requirements Physiological dental chair operated by electricity.

C. Technical Specifications

1. It has double articulated head rest.
2. It has two 3 way syringes (tip autoclavable, with 6 spare tips) one on unit side and other on the assistant side.
3. It has two high speed Air Rotor terminals with water control on coupling.
4. Has one air motor terminal.
5. Arm of unit should be pneumatically locked
6. Removable auxillary tray (stainless steel)
7. It has latest foot operated LED light (min 25,000 LUX)



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395
Email: info@ais-india.com, info@hplcttechnologies.com, info@analyticalgroup.net
Web: www.ais-india.com / www.analyticalgroup.net / www.hplcttechnologies.com

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement

8. It has Rotatable Water System with removable spittoon
9. It has Medium Vacuum Suction
10. It has following programmes –
11. Two programmable working positions
12. Spitting and last working position with light ON and OFF automatically
13. Return to Zero position with light OFF automatically
14. It has option to Lock the movements of chair
15. It has LED based X-ray viewer
16. It is provided with right arm (options for Fixed, Lateral 90 degree swivel available)
17. It is provided with one doctor's stool with adjustable height & backrest tilt including an adjustable ring for foot rest.



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395
Email: info@ais-india.com, info@hplctechnologies.com, info@analyticalgroup.net
Web: www.ais-india.com / www.analyticalgroup.net / www.hplctechnologies.com

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement

Analytical Technologies Limited.

www.analyticalgroup.net