



ADF 2218

## ***Freezer 1 -20°C (12cu.ft) with stabilizers***

Capacity : more than 430 ltrs.

Dimensions : around 600 x 600 x 1300

Insulated inner doors to minimized air loss

Well insulated with 130 mm non CFC, foamed in place polyurethane

Hermetically sealed cooling system with durable compressors

Easy to clean condenser filter

Pressure equalization port to allow easy to access to cabinet after door closer.

Heated door seals (hot gas) to prevent ice build up

Minimum of 15 (2 —) cryo boxes per rack

Minimum of 10 (3 —) cryo boxes per rack

Able to run at 220 Volts and 60 Hz

With safety alarm

All accessories& electric fitting included

Automatic voltage stabilizer be provided