

ANALYTICAL TECHNOLOGIES DEEP FREEZER

Temperature Range :-10°C to -40°C

Special Features:

Upright type (Vertical)

Inside temperature ranging from -10°C to -40°C can be set freely, controlling precision is 1

Stainless steel interior, Outer is painted steel panel

Front opening lockable door with full height handle

110mm foaming insulation 1 set

• Inner 3 shelves Stainless steel with inner door

• imported Danfoss compressor, 2-4 hours cooling temperature come to -40 steady.

• Fan electromotor with internationally famous of Germany EBM

• Microprocessor – based temperature control system, digital temperature display

• Keyboard lock and password protected configuration page

• Audible and visual alarm

High or low temperature alarm,

• Power failure

• Low battery alarm

• Door ajar alarm,

• Filter blocking system failure alarm with temperature printer, record 7 days data.

• 72 hour battery back up for printer and controller display.

• 4 units casters for easy handling

• Mixture gas with R134a, R404a, R23, R24, N2, R600a

• Certificate ISO 9001, ISO 14001, ISO 13485

• Power Supply 220V/50Hz 1 Phase, can change as 220V 60Hz or 110V 50HZ or 110V 60HZ

Optional:

Chart recorder, CO2 back up, storage racks, boxes

Specification:

Capacity : 531 Liters

External Size (W*D*H) MM: 940*840*1999

Inside Size (W*D*H) MM: 650*600*1345

Input Power: 1000W

Weight (Nt/Gt): 235/250kg

