

# Digital Balance, 0.1mg DB-3001 series



## Analytical Balance Electronic Analytical Balance

<b>Capacity</b>	➤	220gm
<b>Readability</b>	➤	0.1mg
<b>Linearity (+/-)</b>	➤	0.2mg
<b>Pan Size</b>	➤	90mm dia
<b>Response Time</b>	➤	2-3 Second
<b>Display</b>	➤	Back Life LCD Display
<b>Calibration</b>	➤	Automatic External
<b>Unit of Measure</b>	➤	G.mg, GN, mo, oz, dwt
<b>Tare Range</b>	➤	Full
<b>Operation Temperature</b>	➤	15 to 40 deg C
<b>Power Supply</b>	➤	AC 230V/115V, 50/60Hz

## Regulatory compliances



## Corporate Social Responsibility



**Analytical  
Foundation**

Analytical Foundation is a nonprofit organization (NGO) founded for the purpose of:

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various fields of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)
2. Improving quality of life by offering YOGA Training courses, Workshops / Seminars etc.
3. ANALYTICAL FOUNDATION aims to DETOXIFY human minds, souls and bodies by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Comps etc.

## Reach us @



**Analytical**®  
Technologies Limited

HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W: [www.ais-india.com](http://www.ais-india.com)  
[www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide