

Digital Conductivity Meter With Cell DCM-3001 Series



Product Description

Table Top Digital Conductivity Meter is used for testing municipal and industrial aqueous solutions. This Conductivity Meter has cell constant adjustment facility.

The Table Top Digital Conductivity Meter has temperature coefficient correction facility. This Conductivity Meter has platinum DIP conductivity cell. This Conductivity Meter has accessories such as dust cover and stand.

The Technical Specifications:

Display: 3 1/2 Digit LED

Temperature Compensation: Manual 0-50 Degree C

Measuring Frequency: 1000 Hz

Power: 230V+/- 10%

Dimensions: 76X 275X 175 mm

Weight: 2 Kg

Regulatory compliances



Corporate Social Responsibility



Analytical
Foundation

Analytical Foundation is a nonprofit organization (NGO) founded for the purpose of:

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various fields of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org
2. Improving quality of life by offering YOGA Training courses, Workshops / Seminars etc.
3. ANALYTICAL FOUNDATION aims to DETOXYFIFY human minds, souls and bodies by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Comps etc.

Reach us @



Analytical[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide