



Dual Dome ceiling suspended operation light (Surgical Light)



ATL 2900.

1. Large dome-100000 lux or more
2. Satellite dome- 95000 lux or more
- B.** Colour (light) temperature- 4200-4500 deg Kelvin
- C.** Type of lamp- Halogen
- D.** Life time of lamp-(1000 to 10000) hours Approx
- E.** Automatic switches to reserve bulb.
- F.** Radiated heat less than 5-8 deg C
- G.** Adjustable position- at least 3 axis up to 60 deg
- H.** Non fogging surface
- I.** Dimming facility of light of 3-6 levels without the change of colour temperature
- J.** Depth of illumination 80-130cm
- K.** Light emission area 1500-2400 cm sq
- L.** Timer indicating number of hours used
- M.** Detachable autoclavable handles - 25 Nos.