



## Analytical Technologies Examination couch



**ATL 4660 ET**

### Specification:-

- A. Overall approx size: 1890 mm L x 560mm W x 840mm H.
- B. Fixed upholstered top 64mm thick in two sections.
- C. Body frame work made from 18G.
- D. CRCA sheet and 20 mm x 40mm x 18 G MS.
- E. Rectangular Tubes Couch fitted with stainless steel Legs.
- F. Headrest adjustable on gas spring.
- G. Upper section of box approx size 1220 mm L x 460 mm W x 630 mm H with three sliding drawers of approx size 320 mm L x 430 mm W x 75 mm H.
- H. Lower section comprises of three cabinets of approx inside size 350 mm L x 440 W mm x 430 H mm with separate doors & lock.
- I. B.P.apparatus tray made of 18 G MS sheet of approx size 350 mm L x 120 mm W X 20 mm H provided on a swinging rod rotating through a bush welded on the body of the couch.

J. Sliding Inbuilt Step Stool.

K. All mild steel components is thoroughly pre-treated chemically to remove rust and foreign matter like Grease, Oil etc. by dip tank process pre-treatment system.

M. The treated Metal Surface has coating of Epoxy Polyester Powder and oven baked at 180 degree to 200 degree

N. Centigrade to avoid contamination of the clean metal surface from dust particles.



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395  
Email: [info@ais-india.com](mailto:info@ais-india.com), [info@hplctechnologies.com](mailto:info@hplctechnologies.com), [info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
Web: [www.ais-india.com](http://www.ais-india.com) / [www.analyticalgroup.net](http://www.analyticalgroup.net) / [www.hplctechnologies.com](http://www.hplctechnologies.com)

**Note :** Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement

**Analytical Technologies Limited.**

**[www.analyticalgroup.net](http://www.analyticalgroup.net)**