

## Analytical Technologies preDIABETIC Glucometer



### SPECIFICATION :

**Measurement Range : 1.1 – 33.3 mmol/L (20-600mg/dL)**

**Blood Sample : Fresh capillary whole blood**

**Test Time: 7 Seconds**

**Operating Temperature : 4 Deg C to 35 Deg C**

**Humidity : < 80%**

**Memory : 220 sets with Date and Time**

**Date Management : 7,14,28 Days Averaging**

**Display:LCD**

**Power Supply : 3V button cell**

**Battery Life : 1)1000test**

2) Replaceable without using any tools

**Size : 85\*55\*17mm**

**Weight : 50g**

**Corporate & Reg. Office:** Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395

Email: [info@ais-india.com](mailto:info@ais-india.com), [info@hplctechnologies.com](mailto:info@hplctechnologies.com), [info@analyticalgroup.net](mailto:info@analyticalgroup.net)

Web: [www.ais-india.com](http://www.ais-india.com) / [www.analyticalgroup.net](http://www.analyticalgroup.net) / [www.hplctechnologies.com](http://www.hplctechnologies.com)

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice.

**We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement.**