



## HOT AIR OVEN

### HAO-2100 Series

#### Specification:

Inner Chamber Size: 455x455x605 mm

Volume: 125 ltrs

Temperature variation (time): +/- °C 0.5

Temperature deviation (spatial): +/- °C 0.5

Readability/ Set ability: °C 0.5

Temperature range: 100-300°C

Degree C above ambient to 250 degree C

Sensor thermocouple: Type J

Control: Nano technology based ergonomic controller

Display: LCD Power : 220 VAC, 1 Phase, 50 Hz

