



Heater plate

Hot Plates have rectangular body made of PCRC sheet finished in enamel paint, rectangular cast iron plate with insulated elements to work on 220v ac 50 hz single phase. Digital temp. Indicator-cum-controller(optional)

Sizes of plates in inches rating:
Size : 12 x 18 x 6" 2.0 kw

Corporate & Regd. Office :

Analytical House, # E67 & E68, Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015, INDIA

Tele : 0265-2253620 / 2252839 Fax : 0265-2254395

E-mail : info@hplctechnologies.com / info@ais-india.com Web : www.ais-india.com / www.analyticalgroup.net