



## Laryngoscope set with neonate blades



ATL 4700 LYS

- A.** Constituted of large hollow, cylindrical, slightly ribbed handle and a set of depressors in Stainless steel.
- B.** Handle is made of either chromium-plated or stainless steel and can be opened at an Extremity to insert two alkaline batteries (size **AA**, 1.5 Volts). The other end has a stud contact which fits the various sizes and types of depressors.
- C.** 2 x straight depressors, Miller type No. **00**, **0**, **1** with halogen bulb
- D.** Presented in suitable protective plastic box.