



Analytical Technologies Monitoring Equipment (ECG, Pulse, Temp, SpO2)



ATL 3501 ME

Specification:-

The monitor should have:

- A. Modular monitor High – resolution colour TFT display of minimum 10" or more.
- B. Able to monitor ECG, NIBP, 2 IBP, SpO2. Temperature and Respiration.
- C. Plethysmograph with perfusion indicator
- D. Monitor seven channel.
- E. 24 Hrs. graphical / tabular trends.
- F. NIBP trends memory is at least 50 readings (tabular).
- G. Suitable for Adult / paediatric/neonate.
- H. Selectable Arrhythmia detection.
- I. Inbuilt three channel recorder.

J. Have Graded and Colour coded alarms

K. User selectable screen formats and user – friendly menu driven functions.

L. Battery backup for at least 3 Hrs.

M. Supplied with: One 3 lead ECG cable, Reusable SpO2 (adult, paediatric, neonate) sensor, NIBP cuffs (each for Adult, child and neonate), IBP cable.



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395
Email: info@ais-india.com, info@hplctechnologies.com, info@analyticalgroup.net
Web: www.ais-india.com / www.analyticalgroup.net / www.hplctechnologies.com

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement

Analytical Technologies Limited.

www.analyticalgroup.net