



MUFFLE FURNANCE

Products | Applications | Software | Accessories | Consumables | Services

MUFFLE FURNANCE

Features :

- Maximum temperature is 1200°C
- œ A choice of temperature controller is available for each size
- œ Double skinned outer case for cool exterior
- œ Automatic switch off when door is opened



Analytical®
Technologies Limited

Analytical Group Companies
Analytical Technologies Limited
Analytical Bio-Technologies
Analytical Bio- Med
HPLC Solutions

Branch Office / Distributor

Corporate & Regd. Office :

Analytical House, # E67 & E68. Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015. INDIA
Tele : 0265-2253620 / 2252839 Fax : 0265-2254395 E-mail : info@hplctechnologies.com / info@ais-india.com
Web : www.ais-india.com / www.analyticalgroup.net