



## Analytical Technologies Nebulizer



**ATL 7401 NBU**

### Specifications:-

- A. Description of Function
  - 1. Nebulizer is a device used to administer medication to people in forms of a liquid mist to the airways. It is commonly used in treating cystic fibrosis, asthma, and other respiratory diseases
- B. Operational Requirements
  - 1. Heavy duty compact Nebuliser.
- C. Technical Specifications
  - 1. Compact, light weight, low noise
  - 2. Durable long life compressor. Suitable for heavy duty/ institutional (hospital) use, is able to run uninterruptedly for one hour, Max Press= 2.0-2.5 bars
  - 3. Produce particle of size 1-5 micron
  - 4. Aluminium cabinet painted with epoxy powder.
  - 5. Piston-type electric aspirator that offers high performance and great durability.
  - 6. Protective thermal cut out relay
  - 7. Air delivery rate app.15 L/min.
  - 8. 24 hours continuous work for hospital use.