



## Analytical Technologies Nursing Trolley



**ATL 4230 NT**

### Specification:-

- A. Size: 750\*450\*900mm (Approx)
- B. Made of high-strength ABS material
- C. One-piece ABS plastic top board with raised-edge design, covered transparent soft plastic glass.
- D. With four plastic-steel columns.
- E. Five ABS drawers: 2 small, 2 middle & 1 big drawers, inner with dividers can be organized easily and freely, label cards for the trolley and drawers.
- F. Stainless steel guard rail
- G. Sliding side shelf (Optional), Writing Shelf(Optional)
- H. Multi bin container, Over bridge tray, Dust basket, Needle disposal holder, IV pole, Utility container, storage box
- I. Centralized lock
- J. Four casters with brakes



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395  
Email: [info@ais-india.com](mailto:info@ais-india.com), [info@hplctechnologies.com](mailto:info@hplctechnologies.com), [info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
Web: [www.ais-india.com](http://www.ais-india.com) / [www.analyticalgroup.net](http://www.analyticalgroup.net) / [www.hplctechnologies.com](http://www.hplctechnologies.com)

**Note :** Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement