

Analytical Technologies Osmo-meter-AO 3600 Series



Specifications

| | |
|--------------------------------|--|
| Model | AO-3600 |
| Cat.No. | 4585 |
| Scale | Urine Osmolality |
| Measurement Range | 0 to 1500 mOsmol/kgH ₂ O 10.0 to 35.0°C |
| Resolution | 10 mOsmol/kgH ₂ O / 0.1°C |
| Measurement Accuracy | ±20 mOsmol/kgH ₂ O (when using alkaline batteries) See "Detailed explanation on Osmolality" / ±1°C |
| Temperature compensation range | 10 to 35°C (Automatic Temperature Compensation) |
| Sample Volume | 0.3ml |
| Measurement Time | 3 seconds |
| Power Supply | Size AAA alkaline Battery x 2 |
| Battery Life | About 11,000 measurements with new alkaline batteries. |
| International Protection Class | IP65 Dust-tight and Protected against water jets. |
| Dimensions & Weight | 55 x31 x109 mm, 100g (Main Unit only) |