



Portable Ultrasound B & W



ATL 3440 PU

Scanning mode: linear and convex
Display mode: B mode, M mode and B+M
Monitor screen be 12" high resolution latest technology
Display depth:25 cm max
Magnification: 8 steps
Facility for image zoom, image reverse, image freeze and image magnification, reduction
Alphanumeric keyboard for character entry and character display
Complete Obstetrics calculation package.
Gynecology measurement facility for uterus, cervix.
probes :
3.5 Mhz convex probe
6 Mhz ICT probe
ultrasound gel -5 Ltr pack with dispensing bottle
Black and white video thermal printer with at least 10 no. Rolls, mobile stand.



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395
Email: info@ais-india.com, info@hplctechnologies.com, info@analyticalgroup.net
Web: www.ais-india.com / www.analyticalgroup.net / www.hplctechnologies.com

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement