



Resuscitation Cart



ATL 5112 RC

Technical Specification:

- A.** Compact unit with minimum 'foot print' and provided with good quality castors with braking system.
- B.** Modular storage system, rust proof, allows full extension drawers.
- C.** 40-45" height
- D.** Swing outside pod left side for extra space with 4 Nos. tilt out bins & one pod tray for storage.
- E.** Short utility hooks to hang items, with one utility pole and clamps.
- F.** Single over bridge with minimum 2 hangers rails to fit assorted baskets, shelves & bins.
- G.** Specimen bag holder along with over bridge with clamps.
- H.** Utility bin on lower hanger rail and 2 utility bins on upper hanger rail with three tape/label dispenser.
- I.** Multiple drawers 3" drawers 9" drawers which are interchangeable with central locking system.
- J.** The firm clearly indicates in the technical bid itself that the prices of all standard accessories and the accessories mentioned in the specification are included in the quoted price.



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395
Email: info@ais-india.com, info@hplctechnologies.com, info@analyticalgroup.net
Web: www.ais-india.com / www.analyticalgroup.net / www.hplctechnologies.com

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement