



Resuscitation Equipment/bags (Adult)



ATL 4800 SERIES

A. To have Retromolar Intubation fibroscope for unexpected difficult airways.

1. Tip Distal Bending 40°.
2. To be movable eyepiece
3. To have a light source connection
4. With length 40-42cms and dia 5-6 mms.
5. ET tube holder be provided
6. take min. 5.5 size of ET tube

B. Portable LED light source be provided with;

1. Illumination not less than 50000 Lux
2. run on two 3V photo batteries
3. Burning life be more than 100 minutes
4. Ergonomically designed and can be connected to both the fibrescopes
5. Life of LED be close to 50000 hrs

C. One Laryngoscope with rechargeable battery pack and blade with fiberoptic mechanism be provided to be used on both adult and pediatric patients with charger.

D. Other accessories like, magill forceps be provided.



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395
Email: info@ais-india.com, info@hplctechnologies.com, info@analyticalgroup.net
Web: www.ais-india.com / www.analyticalgroup.net / www.hplctechnologies.com

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement



- E.** have Emergency Cricothyroidotomy for pediatric and adult
 - 1. Disposable blades
 - 2. Dialator
- F.** has Combitube size 37Fr. With complete kit

- G.** have Intubating Laryngeal Mask Airways with Following Components:
 - 1. ILMA Sizes 3 & 4.
 - 2. ILMA Tubes ID 7mm & 7.5mm.
 - 3. Tube Stabilizing rod
 - 4. Cuff deflator
- H.** have Laryngeal Mask Airways Sizes 1, 2 and 4
- I.** Handy and strong brief case/bag be provided to keep all the instruments safe.
- J.** Set of disposable percutaneous tracheotomy kit for adult and pediatric.
- K.** have standard AMBU bag for pediatric and adult.
- L.** Mechanical suction pump with suction catheter and stomach tubes.
- M.** have Aluminum Oxygen reservoir 2 Liter with oxygen tube and cather.
- N.** Oxygen pressure reducer, regulable 0-15 liter with coupler for respirator.
- O.** Ventilating bag
- P.** Lubricant
- Q.** Blood pressure meter, boso K-II
- R.** Stethoscope
- S.** Rescue blanket gold/silver
- T.** Infusion system.