



**Analytical**<sup>®</sup>  
Technologies Limited



### SOXHLET EXTRACTION SET (ATL-SEU)

Standard mild steel body and upright fittings. Body finished in stove heating paint and other parts are duly chrome plated. Housing provided with temp. up to 300 °C . Provided with a common energy regulator to control all units, individual on-off switches & pilot lamps. To work on 220VAC

The unit of 6 mantels rating 750 watts.heating mantels for surface



**Analytical**<sup>®</sup>  
Technologies Limited

Analytical Group Companies  
Analytical Technologies Limited  
Analytical Bio-Technologies  
Analytical Bio- Med  
HPLC Solutions

Branch Office / Distributor

Corporate & Regd. Office :

Analytical House, # E67 & E68. Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015. INDIA  
Tele : 0265-2253620 / 2252839 Fax : 0265-2254395 E-mail : info@hplctechnologies.com / info@ais-india.com  
Web : www.ais-india.com / www.analyticalgroup.net