



Stretcher Trolley



ATL 5600 SW

- A.** 213L x 56W x 81H cms.
- B.** Frame work made of 31 mm x 1.60 mm vertical & 25 mm x 1.22 mm horizontal CRC Tubes Trolley mounted on 15 cms dia castors – 2 with brakes.
- C.** Removable stretcher top made of 1.22 mm aluminium sheet with S.S. handle at both end with 25mm thick suitable rubber.
- D.** Mattress covered with good quality rexine.
- E.** All mild steel components are thoroughly pre-treated chemically to remove rust and foreign matter like Grease; Oil etc. by diptank process pre-treatment System.
- F.** The treated Metal Surface have coating of Epoxy Polyester Powder and oven baked at 180 degree to 200 degree Centigrade to avoid contamination of the clean metal surface from dust particles.