



VACUUM EXTRACTOR



ATL 4420 VE

- A.** Noiseless suction unit have fast vacuum built up
- B.** Vacuum have maximum -90 k pa/-675 mm Hg, Suction capacity 50ltr/min/1kg/cm² & bottle capacity 3ltrs
- C.** Suction system have piston / cylinder (self lubricating)
- D.** have mechanical overflow protection system.
- E.** Set of silicon cups 50mm & 60mm-2each
- F.** Set of bird cups, stainless steel 40mm, 50mm & 60mm-2each
- G.** Machine can be operated on 220-240 V AC single phase
- H.** be provided with a gauge to display of vacuum generated
- I.** have vacuum release valve