



**Analytical**<sup>®</sup>  
Technologies Limited



## VORTEX SHAKER

### SPECIFICATION

Useful for accelerated mixing of solutions in test tubes, small flasks or centrifugal tubes etc. The housing is 150x150 mm The speed is regulated by a speed control devices. Shorts eccentric stroke motor imparts vigorous agitation to the tubes or flasks . Body is strong and finished in white stoving enamel/pow der coated paint with mat finished colour combinations.The equipment is workable on 220VAC 50 Hz single phase.



**Analytical**<sup>®</sup>  
Technologies Limited

Analytical Group Companies  
Analytical Technologies Limited  
Analytical Bio-Technologies  
Analytical Bio- Med  
HPLC Solutions

Branch Office / Distributor

Corporate & Regd. Office :

Analytical House, # E67 & E68. Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015. INDIA  
Tele : 0265-2253620 / 2252839 Fax : 0265-2254395 E-mail : info@hplctechnologies.com / info@ais-india.com  
Web : www.ais-india.com / www.analyticalgroup.net