



Analytical[®]
Technologies Limited



Vicat Needle Apparatus

Specification

The apparatus consists of a metallic frame bearing a freely movable rod with a cap at top, one vicat mould. Split type and glass base part and one set of needles one each initial needle and consistency plunger



Analytical[®]
Technologies Limited

Analytical Group Companies
Analytical Technologies Limited
Analytical Bio-Technologies
Analytical Bio-Med
HPLC Solutions

Branch Office / Distributor

Corporate & Regd. Office :

Analytical House, # E67 & E68. Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015. INDIA
Tele : 0265-2253620 / 2252839 Fax : 0265-2254395 E-mail : info@hplctechnologies.com / info@ais-india.com
Web : www.ais-india.com / www.analyticalgroup.net