



Video Gastroscope



GASTRO 3291V SERIES

A. Other Accessories:

1. Slim, light weight, fully immersible in disinfectant solution
2. Field of view: 120° to 150°
3. Depth of field: 3 - 100 mm
4. Tip Deflection: Up/Down 210 °/90°
5. Left/Right 100°/100°
6. Rigid distal diameter: Less than 9.8 mm
7. Insertion tube diameter: Not more than 9.8 mm
8. Instrument channel: 2.8 to 3.2 mm
9. Working length: 1050 to 1100 mm
10. Total length: 1350 to 1400 mm
11. Xenon light source cum video processor
12. Integrated or separate units with Xenon lamp
13. Color System : Single CCD Color
14. Lamp: 300 W Xenon
15. Appropriate connectors
16. Image display size: Full and small screen display on monitor
17. Video outputs: RGB, Y/C & composite output
18. Light control system: Automatic and manual control



Corporate & Reg. Office: Analytical House, E67 & 68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda, Gujarat-390 015, INDIA. Tel.: 0091- 265-2253620, 2252839 Telefax: 0091-265-2254395
Email: info@ais-india.com, info@hplctechnologies.com, info@analyticalgroup.net
Web: www.ais-india.com / www.analyticalgroup.net / www.hplctechnologies.com

Note : Company reserves right to add/delete/modify the content / technical specification of the catalogue without prior notice. We are the only Tailor Making Company who provide Analytical Equipment based on Customer Requirement



19. Cooling system: Forced air cooling
20. Video Processing System: (common for all videoscope)
21. The Video Processor have facility of compact flash memory card. Processor have facility to connect balloon enteroscope and endoscopy ultra sound.
22. LCD Digital Monitor 21", (Colour medical grade)
23. Endoscopy trolley (from OEM): 1
24. Software for image capturing. DVD recorder, color printer, Archiving memory, report generator.
25. PC – Latest configuration with HDD 500 GB, RAM 4GB, TFT 21 inch, DVD R/W Drive

B. Other Accessories:

1. Biopsy forceps
2. Bite block
3. Leakage tester
4. Cleaning brush
5. Aspiration needle
6. EB Forceps