

Water Distillatory WD 3000 Series



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

Product Description

Uses:

The series of apparatus have tap water as source to produce pure water by electrical heating distilling. It is applied in health and medicine units, chemical industries, scientific research institutions and labs etc.

Characteristics:

1. Made of high-quality stainless steel by stamping and welding.
2. Characterized by anti-corrosion, age-resistant, easy operation and stable function, and safety and durability.
3. Coiled stainless steel tube condenser with good heating exchange and large water output.
4. Special water level design, under the low water level condition, the alarm system will work and cut off the power supply quickly. This makes sure the heating element no damage.
5. Automatic water supply function, when the water level is low, the floater will automatically decline, the water comes into the equipment ensure continue working, save time and make sure the high safety.
6. The 5 Liters Water Distiller is equipped with the standard Type.
The 10 Liters and 20 Liters Water Distiller are standard Type.
7. Each water distiller Stainless Steel Tank to storage the distilled water. The possibility is connection to system of the centralized water supply.
The 5Liters Water Distiller has one 10 Liters Stainless Steel Tank, both the 10 Liters and 20 Liters Water Distiller have one 20 Liters Stainless Steel Tank.8.. Warranty 24 months from the date of putting into operation.
8. ISO9001:2008; CE; ISO 13485:2003.

Supply voltage	Output L/h	Heating power (KW)	Rough weight (kg)	Net weight (kg)	Overall product dimension (mm)	Overall package dimension (mm)
380	10	7.5	7.5	10.5	350×270×830	390×350×890

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

- HPLC Servicing** :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.
- Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.
- AMC's/CMC** :AMC's/CMC :We offer user training both in-House and at customer sites on PLC principles, operations, troubleshooting.
- Validations** :Validations :We have protocols for carrying out periodic Validations as per GLP/ GMP/U SFDA norms.
- Instruments** :Instruments :We offer instruments / Renting Services Modules like pumps, detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies

						
UV/VIS Spectro 2060+ Double Beam	FTIR	Gas Chromatograph 3000	Gas Chromatograph 2979 Plus	Flash Chromatograph	Atomic Absorption Spectrophotometer	Liquid Partial Counter
						
Optical Emission Spectrophotometer	DSC/TGA	NOVA 2020 plus Automated Bio Chemistry Analyzer	HEMA 2020 Hematology Analyzer	Micro Plate Reader/Washer	URINOVA 2800 Urine Analyzer	Total Organic Carbon
						
Fully Automated CLIA	NOVA Basic Semi-Auto Chemistry Analyzer	PCR/Gradient PCR/ RTPCR	Blood Gas Analyzer	Random access Analyzer for immunoassay Proteins & clinical chemistry	Semen Analyzer	Water purification system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.