



Analytical[®]
Technologies Limited



Laboratory Refrigerator

SPECIFICATION

Temperature Sensitivity: \pm °C 1 or better

Spatial Deviation In Temperature: \pm °C 1 or better

Readability °C : 0.1

Temperature range °C: 4°C to 8°C

Temperature Controller :Solid State digital Controller/PID.

Display :LED/LCD

Capacity: 15 CFT



Analytical®
Technologies Limited

Analytical Group Companies
Analytical Technologies Limited
Analytical Bio-Technologies
Analytical Bio- Med
HPLC Solutions

Branch Office / Distributor

Corporate & Regd. Office :

Analytical House, # E67 & E68. Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015. INDIA
Tele : 0265-2253620 / 2252839 Fax : 0265-2254395 E-mail : info@hplctechnologies.com / info@ais-india.com
Web : www.ais-india.com / www.analyticalgroup.net